



HANDY
MASSAGER

Rilamingo

Instruction Manual

For
storage

Medical device: home electric massager


Lourdes Handy Massager "Flamingo"
AX-KXL3200



*Please understand some designs and specifications may have changed.


Thank you very much for purchasing AX-KXL3200. Before use, read this instruction manual and use the product accordingly. Please keep this manual with care after reading it.

Name of the parts




Preparation before use


How to put batteries



STEP 4 Put battery case to main body.



STEP 5 Lock: Align the ● mark on the battery cap with the lock mark as figure.



How to use

STEP 1 Put vibration part to body. Start massaging.

example

[shoulder]



[lower back]



[palm]



STEP 2 Release vibration part from body. Stop massaging.

Compatible with waterproof function, IPX5.

*IPX5: There is no effect on water jets from all directions.



Prohibited

*Do not soak, use in bathroom or wash with water.
*Do not use to thin part of the skin like head, knees or elbows, weak part of skin, face, throat and belly.
< It may cause an accident or injury.>